

Stoner Python Package

C.S. Allen, M. Newman, R. Temple, S. Morely and G. Burnell

November 16, 2012

Contents

1	Introduction	3
1.1	Getting the Stoner Package	3
1.2	Using the Stoner Package	3
1.3	Documentation	3
2	Users' Guide	4
2.1	Loading a data file	4
2.1.1	Loading Data from a strong or iterable object	6
2.2	Examining Some Data	6
2.2.1	Data, Column headers and metadata	6
2.2.2	Masked Data and Why You Care	6
2.2.3	Working with columns of data	7
2.2.4	Working with complete rows of data	7
2.2.5	Manipulating the metadata	8
2.2.6	More on Indexing the data	8
2.2.7	Selecting Individual rows and columns of data	9
2.2.8	Find out more about the data	10
2.2.9	Copying Data	10
2.3	Modifying Data	10
2.3.1	Appending data	10
2.3.2	Working with Columns of Data	11
2.3.3	Rearranging Columns of Data	11
2.3.4	Renaming Columns of Data	11
2.3.5	Inserting Columns of Data	12
2.3.6	Deleting Rows of Data	12
2.3.7	Deleting Columns of Data	12
2.3.8	Sorting Data	12
2.4	Saving Data	13
2.5	Plotting Data	13
2.5.1	Plotting 2D data	13
2.5.2	Plotting 3D Data	14
2.5.3	Getting More Control on the Figure	15
2.6	Manipulating and Analysing Data	15
2.7	Manipulating Data	16
2.8	Curve Fitting	17
2.8.1	Simple polynomial Fits	17
2.8.2	Simple function fitting	17

2.8.3	Fitting with limits	17
2.9	More AnalyseFile Functions	19
2.9.1	Basic Data Inspection	19
2.9.2	Thresholding and Interpolating Data	19
2.9.3	Smoothing and Differentiating Data	20
2.9.4	Peak Finding	20
2.10	Working with Lots of Files	20
2.10.1	Getting a List of Files	20
2.10.2	Doing Something With Each File	21
2.10.3	Sorting, Filtering and Grouping Data Files	21
3	Cookbook	23
3.1	Extract X-Y(Z) from X-Y-Z data	23
3.2	Mapping X-Y-Z data to Z(X,Y) data	23
4	Developer's Guide	24
4.1	Adding New Data File Types	24

1 Introduction

This manual provides a user guide and reference for the Stoner python package. The Stoner python package provides a set of python classes and functions for reading, manipulating and plotting data acquired with the lab equipment in the Condensed Matter Physics Group at the University of Leeds.

1.1 Getting the Stoner Package

The source code for the Stoner python module is kept in CVS revision control on the stonerlab server. A stable release of the code is available for copying and use in `\\stonerlab\data\software\python\stable\`. The development code can be obtained by checking out the PythonCode module with a CVSROOT of `:ext:cvs@stonerlab.leeds.ac.uk:/home/cvs/`. Appropriate ssh keys for the cvs user account are kept in `\\stonerlab\data\software\CVS\`.

The Stoner Package currently depends on a number of other modules. These are installed on the lab machines that have Python installed. Primarily these are Numpy, SciPy and Matplotlib. Windows installable versions are kept in `\\stonerlab\data\software\Python\`. The easiest way to get a Python installation with all the necessary dependencies for the Stoner Package is to install the *Enthought Python Distribution*. Windows install file are kept in `\\stonerlab\data\software\python`

1.2 Using the Stoner Package

The easiest way to use the Stoner Package is to add the path to the directory containing Stoner.py to your PYTHONPATH environment variable. This can be done on Macs and Linux by doing:

```
cd <path to PythonCode directory>/src
export PYTHONPATH='pwd':$PYTHONPATH
```

On a windows machine the easiest way is to create a permanent entry to the folder in the system environment variables. Go to Control Panel -> System -> Advanced Tab -> click on Environment button and then add or edit an entry to the system variable PYTHONPATH.

Once this has been done, the Stoner module may be loaded from python command line:

```
>>> import Stoner

or

>>> from Stoner import *
```

1.3 Documentation

This document provides a user guide to the Stoner package and its various modules and classes. It is not a reference to the library but instead aims to explain the various operations that are possible and provide short examples of use. For the API reference for the library, please see the *Python Code API* compiled windows help file. Rowan has also made a single sided “cheat sheet” that summarises the examples in this user guide.



The code is still under active development to fix bugs and add features. Generally things don't get deliberately broken, but accidents happen, so if

something stops working, please either fix and commit the code or tell Gavin.

2 Users' Guide

The Users' Guide provides a brief overview of the functions contained within the Stoner module and so basic examples of how the module can be used.

The Stoner module provides several Python classes that can be used to manipulate experimental data. The main class that provides the basic functionality is the `DataFile` class. This handles loading data, finding and manipulating meta data, selecting rows or columns of data, adding or removing data, and saving data.

The `PlotFile` class is a descendent of `DataFile`, meaning it shares all the same functionality as `DataFile`, but in addition has methods to present data graphically. The `AnalyseFile` class is another descendent of `DataFile`, but provides extra methods to fit curves, smooth and differentiate data, find peaks and carry out other simple analysis operations.

2.1 Loading a data file

The first step in using the Stoner module is to load some data from a measurement.

```
>>>import Stoner
>>>d=Stoner.DataFile('my_data.txt')
>>>d=Stoner.VSMFile('my_VSM_data.fld')
```

In this example we have loaded data from `my_data.txt` which should be in the current directory – here we are assuming that `my_data.txt` contains data in the *TDI Format 1.5* which is produced by the LabVIEW rigs. Assuming that the file successfully loads, `d`, is an instance of the `DataFile` object. Here the `DataFile` constructor has been used to both create the instance and load the data in one go.

The second example shows the use of one of the sub-classes of the `DataFile` object to load data from a specific instrument (in this case the VSM).

This is an API change from earlier versions of the Stoner package where a second parameter on the constructor of the `DataFile` object was used to identify the type of data file. This syntax is now depreciated !).



The possible sub-classes are:

DataFile Tagged Data Interchange Format 1.5 – the default format produced by the LabVIEW measurement rigs

VSMFile The text files produced by the group's Oxford Instruments VSM

BigBlueFile Datafiles produced by VB Code running on Big Blue. The *BigBlue* version of the `DataFile.load` and `DataFile` constructors takes two additional parameters that specify the row on which the column headers will be found and the row on which the data starts.

CSVFile Reads a generic comma separated value file. The **CSVFile** load routine takes four additional parameters to the constructor and load methods. In addition to the two extra arguments used for the *BigBlue* variant, a further two parameters

specify the delimiters for the data and header rows. **CSVFile** also offers a **save** method to allow data to be saved in a simple delimited text way (see Section 2.4 for details).

XRDFFile Loads a scan file produced by Arkengarthdale - the group's Brucker XRD Machine.

SPCFile Loads a Raman scan file (.spc format) produced by the Renshaw and Horiba Raman spectrometers. This may also work for other instruments that produce spc files, but has not been extensively tested.

BNLFile Loads a SPEC file from Brookhaven (so far only tested on u4b files but may well work with other synchrotron data). Produces metadata Snumber: Scan number, Stype: Type of scan, Sdatetime: date time stamp for the measurement, Smotor: z motor position.

TDMSFile Loads a file saved in the National Instruments TDMS format

QDSquidVSMFile Loads data from a Quantum Design SQUID VSM as used on the I10 Beamline in Diamond.

OpenGDAFile Reads a scan file generated by OpenGDA – a software suite used for synchrotrons such as Diamond.

RasorFile Simply an alias for OpenGDAFile used for the RASOR instrument on I10 at Diamond.

FmokeFile Loads a file from Dan Allwood's Focussed MOKE System in Sheffield.

```
>>>import Stoner
>>>d=Stoner.DataFile()
>>>d.load('my_data.txt')
>>> v=Stoner.VSMFile()
>>>v.load('my_VSM_data.fld')
>>> c=Stoner.CSVFile()
>>> c.load('data.csv',1,0,',',',',',')
```



*The load method, like many of the DataFile methods returns a copy of the Datafile object **as well as** modifying the object itself. The advantage of this is that it is then possible to chain several methods into one command*

Sometimes you won't know exactly which subclass of **DataFile** is the one to use. Unfortunately, there is no sure fire way of telling, but **DataFile.load** will try to do the best it can and will try all of the subclasses in memory in turn to see if one will load the file without throwing an error. If this succeeds then the actual type of file that worked is stored in the metadata of the loaded file.



The automatic loading assumes that each load routine does sufficient sanity checking that it will throw an error if it gets bad data. Whilst one might wish this was always true it relies on whoever writes the load method to make sure of this ! If you want to stop the automatic guessing from happening use the **auto_load=False** keyword in the load method.

2.1.1 Loading Data from a strong or iterable object

In some circumstances you may have a string representation of a **DataFile** object and want to transform this into a proper **DataFile** object. This might be, for example, from transmitting the data over a network connection or receiving it from another program. In these situations the *left shift operator*, `<<`, can be used.

```
>>> data=Stoner.Core.DataFile() << string_of_data
>>> data=Stoner.Core.DataFile() << iterable_object
```

The second example would allow any object that can be iterated (i.e. has a *next()* method that returns lines of the data file, to be used as the source of the data. The **Stoner.Core.DataFile()** creates an empty object so that the left shift operator calls the method in **DataFile** to read the data in. It also determines the type of the object **data**. This also provides an alternative syntax for reading a file from disk:

```
>>> data=Stoner.Core.DataFile()<<open("File on Disk.txt")
```

2.2 Examining Some Data

2.2.1 Data, Column headers and metadata

Having loaded some data, the next stage might be to take a look at it. Internally, data is represented as a 2D numpy masked array of floating point numbers, along with a list of column headers and a dictionary that keeps the metadata and also keeps track of the expected type of the metadata (*i.e.* the meta-metadata). These can be accessed like so:

```
>>>d.data
>>>d.column_headers
>>>d.metadata
```

2.2.2 Masked Data and Why You Care

Masked data arrays differ from normal data arrays in that they include an option to mask or hide individual data elements. This can be useful to temporarily discount parts of your data when, for example, fitting a curve or calculating a mean value or plotting some data. One could, of course, simply ignore the masking option and use the data as is, however, masking does have a number of practical uses.

The data mask can be accessed via the *mask* attribute of **DataFile**:

```
>>> import numpy.ma as ma
>>> print d.mask
>>> d.mask=False
>>> d.mask=ma.nomask
>>> d.mask=numpy.array([[True, True, Fale,...False],...,
                        [False,True,...True]])
>>> d.mask=lambda x: x[0]<50
>>> d.mask=lambda x:[y<50 for y in x]
```

The first line is simply the import statement for the numpy masked arrays in order to get the *nomask* symbol. The second line will simply print the current mask. The next two examples will unmask all the data *i.e.* make the values visible and useable. The next example illustrates using a numpy array of booleans to set the mask - every element in the mask array that evaluates as a boolean True will be masked and every

False value unmasked. So far the semantics here are the same as if one had accessed the mask directly on the data via `d.data.mask` but the final two examples illustrate an extension that setting the **DataFile** mask attribute allows. If you pass a callable object to the mask attribute it will be executed, passing each row of the data array to the user supplied function as a numpy array. The user supplied function can then either return a single boolean value – in which case it will be used to mask the entire row – or a list of booleans to mask individual cells in the current row.

By default when the **DataFile** object is printed or saved, data values that have been masked are replaced with a “fill” value of 10^{20} .



This is somewhat dangerous behaviour. Be very careful to remove a mask before saving data if there is any chance that you will need the masked data values again later !

2.2.3 Working with columns of data

This is all very well, but often you want to examine a particular column of data or a particular row:

```
>>>d.column(0)
>>>d.column('Temperature')
>>>d.column(['Temperature',0])
```

In the first example, the first column of numeric data will be returned. In the second example, the column headers will first be checked for one labeled exactly *Temperature* and then if no column is found, the column headers will be searched using *Temperature* as a regular expression. This would then match *Temperature (K)* or *Sample Temperature*. The third example results in a 2 dimensional numpy array containing two columns in the order that they appear in the list (*i.e.* not the order that they are in the data file). For completeness, the **DataFile.column** method also allows one to pass slices to select columns and should do the expected thing.

There is a convenient shortcut for working with cases where the column headers are not the same as the names of any of the attributes of the **DataFile** object:

```
>>>d.Temperature
>>>d.column('Temperature')
```

both return the same data.

2.2.4 Working with complete rows of data

Rows don't have labels, so are accessed directly by number:

```
>>>d[1]
>>>d[1:4]
```

The second example uses a slice to pull out more than one row. This syntax also supports the full slice syntax which allows one to, for example, decimate the rows, or directly pull out the last few rows in the file.

2.2.5 Manipulating the metadata

What happens if you use a string and not a number in the above examples ?

```
>>>d['User']
```

in this case, it is assumed that you meant the metadata with key *User*. To get a list of possible keys in the metadata, you can do:

```
>>>d.dir()
>>>d.dir('Option\:.*)
```

In the first case, all of the keys will be returned in a list. In the second, only keys matching the pattern will be returned – all keys containing *Option:.* For compatibility with normal python semantics: `>>>d.keys()` is synonymous with `>>>d.dir()`.

We mentioned above that the metadata also keeps a note of the expected type of the data. You can get at the metadata type for a particular key like this:

```
>>>d.metadata.type('User')
```

to get a dictionary of all of the types associated with each key you could do:

```
>>>dict(zip(d.dir(),d.metadata.type(d.dir())))
```

but an easier way would be to use the **typeHintedDict.types** attribute:

```
>>>d.metadata.types
```

2.2.6 More on Indexing the data

There are a number of other forms of indexing supported for **DataFile** objects.

```
>>>d[10,0]
>>>d[0:10,0]
>>>d[10,'Temp']
>>>d[0:10,['Voltage','Temp']
```

The first variant just returns the data in the 11th row, first column (remember indexing starts at 0). The second variant returns the first 10 values in the first column. The third variant demonstrates that columns can be indexed by string as well as number, and the last variant demonstrates indexing multiple rows and columns – in this case the first 10 values of the Voltage and Temp columns.

You might think of the data as being a list of records, where each column is a field in the record. Numpy supports this type of structured record view of data and the **DataFile** object provides the *DataFile.records* attribute to do this. This read-only attribute is just providing an alternative view of the same data.

```
>>> d.records
```


2.2.7 Selecting Individual rows and columns of data

Many of the function in the Stoner module index columns by searching the column headings. If one wishes to find the numeric index of a column then the **DataFile.find_col** method can be used:

```
>>> index=d.find_col(1)
>>> index=d.find_col('Temperature')
>>> index=d.find_col('Temp.*')
>>> index=d.find_col(1:10:2)
>>> index=d.find_col(['Temperature',2,'Resistance'])
```

DataFile.find_col takes a number of different forms. If the argument is an integer then it returns (trivially) the same integer, a string argument is first checked to see if it exactly matches one of the column headers in which case the number of the matching column heading is returned. If no exact match is found then a regular expression search is carried out on the column headings. In both cases, only the first match is returned. The final two examples given above both return a list of indices, firstly using a slice construct - in this case the result is trivially the same as the slice itself, and in the last example by passing a list of column headers to look for.

This is the function that is used internally by **DataFile.column**, **DataFile.search** etc and for this reason the trivial integer and slice forms are implemented to allow these other functions to work with multiple columns.

Sometimes you may want to iterate over all of the rows or columns in a data set. This can be done quite easily:

```
>>>for row in d.rows():
.....print row
.....
>>>for column in d.columns():
.....print column
.....
```

The first example could also have been written more compactly as:

```
>>>for row in d:
.....print row
.....
```

In many cases you do not know which rows in the data file are of interest - in this case you want to search the data.

```
>>>d.search('Temperature',4.2)
>>>d.search('Temperature',4.2,['Temperature','Resistance'])
>>>d.search('Temperature',lambda x,y: x>10 and x<100)
>>>d.search('Temperature',lambda x,y: x>10 and
          x<1000 and y[1]<1000,['Temperature','Resistance'])
```

The general form is

```
DataFile.search(<search column>,<search term>[,<list of return columns>])
```

The first example will return all the rows where the value of the *Temperature* column is 4.2. The second example is the same, but only returns the values from the *Temperature*, and *Resistance* columns. The rules for selecting the columns are the same as for

the `DataFile.column` method above – strings are matched against column headers and integers select column by number.

The third and fourth examples above demonstrate the use of a function as the search value. This allows quite complex search criteria to be used. The function passed to the search routine should take two parameters – a floating point number and a numpy array of floating point numbers and should return either *ture* or *False*. The function is evaluated for each row in the data file and is passed the value corresponding to the search column as the first parameter while the second parameter contains a list of all of the values in the row to be returned. If the search function returns `True`, then the row is returned, otherwise it isn't. In thr last example, the final parameter can either be a list of columns or a single column. The rules for indexing columns are the same as used for the `DataFile.find_col` method.

2.2.8 Find out more about the data

Another question you might want to ask is, what are all the unique values of data in a given column (or set of columns). The Python numpy package has a function to do this and we have a direct pass through from the `DataFile` object for this:

```
>>> d.unique('Temp')
>>> d.unique(column,return_index=False, return_inverse=False)
```

The two optional keywords cause the numpy routine to return the indices of the unique and all non-unique values in the array. The column is specified in the same way as the `DataFile.column` method does.

2.2.9 Copying Data

One of the characterisitics of Python that can confuse those used to other programming languages is that assignments and argument passing is by reference and not by value. This can lead to unexcted results as you can end up modifying variables you were not expecting ! To help with creating genuine copies of data Python provides the `copy` module. Whilst this works with `DataFile` objects, for convenience, the `DataFile.clone` attribute is provided to make a deep copy of a `DataFile` object.

This is an attribute not a method, so there are no brackets here !



```
>>> t=d.clone
```

2.3 Modifying Data

2.3.1 Appending data

The simplest way to modify some data might be to append some columns or rows. The Stoner mpmodule redefines two standard operators, `+` and `&` to have special meanings:

```
>>>a=Stoner.DataFile('some_new_data.txt')
>>>add_rows=d+a
>>>add_columns=d&a
```

In these example, *a* is a second DataFile object that contains some data. In the first example, a new DataFile object is created where the contents of *a* are added as new rows after the data in *d*. Any metadata that is in *a* and not in *d* are added to the metadata as well. There is a requirement, however, that the column headers of *d* and *a* are the same – *i.e.* that the two DataFile objects appear to represent similar data.

In the second example, the data in *a* is added as new columns after the data from *d*. In this case, there is a requirement that the two DataFile objects have the same number of rows.

These operators are not limited just to DataFile objects, you can also add numpy arrays to the DataFile object to append additional data.

```
>>>import numpy as np
>>>x=np.array([1,2,3])
>>>new_data=d+x
>>>y=np.array([1,2,3],[11,12,13],[21,22,23],[31,32,33])
>>>new_data=d+y
>>>column=d.column[0]
>>>new_data=d&column
```

In the first example above, we add a single row of data to *d*. This assumes that the number of elements in the array matches the number of columns in the data file. The second example is similar but this time appends a 2 dimensional numpy array to the data. The third example appends a numpy array as a column to *d*. In this case the requirement is that the numpy array has the same or fewer rows of data as *d*.

2.3.2 Working with Columns of Data

2.3.3 Rearranging Columns of Data

Sometimes it is useful to rearrange columns of data. **DataFile** offers a couple of methods to help with this.

```
>>> d.swap_column(('Resistance','Temperature'))
>>> d.swap_column(('Resistance','Temperature'),headers_too=False)
>>> d.swap_column([(0,1),('Temp','Volt')],(2,'Curr'))
>>> d.reorder([1,3,'Volt','Temp'])
>>> d.reorder([1,3,'Volt','Temp'],header_too=False)
```

The **swap** method takes either a tuple of column names/indices or a list of such tuples and swaps the columns accordingly, whilst the **reorder** method takes a list of column labels/indices and constructs a new data matrix out of those columns in the new order. The *headers_too=False* options, as the name suggests, cause the column headers not be rearranged.

2.3.4 Renaming Columns of Data

As a convenience, **DataFile** also offers a useful method to rename data columns:

```
>>> d.rename('old_name','new_name')
>>> d.rename(0,'new_name')
```

Alternatively, of course, one could just edit the `column_headers` attribute.

2.3.5 Inserting Columns of Data

The append columns operator `&` will only add columns to the end of a dataset. If you want to add a column of data in the middle of the data set then you should use the **add_column** method.

```
>>>d.add_column(numpy.array(range(100)), 'Column Header')
>>>d.add_column(numpy.array(range(100)), 'Column Header', Index)
>>>d.add_column(lambda x: x[0]-x[1], 'Column Header', func_args=None)
```

The first example simply adds a column of data to the end of the dataset and sets the new column headers. The second variant inserts the new column before column *Index*. *Index* follows the same rules as for the **DataFile.columnn()** method. In the third example, the new column data is generated by applying the specified function. The function is passed a single row as a 1D numpy array and any of the keyword, argument pairs passed in a dictionary to the optional *func_args* argument.

The **DataFile.add_column** method returns a copy of the DataFile object itself as well as modifying the object. This is to allow the method to be chained up with other methods for more compact code writing.

2.3.6 Deleting Rows of Data

Removing complete rows of data is achieved using the **DataFile.del_row** method.

```
>>>d.del_rows(10)
>>>d.del_rows('X Col', value)
>>>d.del_rows('X Col', lambda x,y:x>300)
```

The first variant will delete row 10 from the data set (where the first row will be row 0). You can also supply a list or slice to **DataFile.del_rows** to delete multiple rows.

If you do not know in advance which row to delete, then the second and third variants provide more advanced options. The second variant searches for and deletes all rows in which the specified column contains *value*. The third variant selects which rows to delete by calling a user supplied function for each row. The user supplied function is the same in form and definition as that used for the **DataFile.search** method.

2.3.7 Deleting Columns of Data

Deleting whole columns of data can be done by referring to a column by index or column header - the indexing rules are the same as used for the **DataFile.column** method.

```
>>>d.del_column('Temperature')
>>>d.del_column(1)
```

2.3.8 Sorting Data

Data can be sorted by one or more columns, specifying the columns as a number or string for single columns or a list or tuple of strings or numbers for multiple columns. Currently only ascending sorts are supported.

```
>>>d.sort('Temp')
>>>d.sort(['Temp', 'Gate'])
```

2.4 Saving Data

Only saving data in the *TDI* format and as comma or tab delimited formats is supported.



The **CSVFile** comma or tab delimited files discard all metadata about the measurement. You absolutely must not use this as your primary data format – always keep the *TDI* format files as well.

```
>>>d.save()
>>>d.save(filename)
>>>d=Stoner.CSVFile(d)
>>>d.save()
>>>d.save(filename,'\t')
```

In the first case, the filename used to save the data is determined from the filename attribute of the **DataFile** object. This will have been set when the file was loaded from disc.

If the filename attribute has not been set *e.g.* if the **DataFile** object was created from scratch, then the **DataFile.save** method will cause a dialog box to be raised so that the user can supply a filename.

In the second variant, the supplied filename is used and the filename attribute is changed to match this *i.e.* **d.filename** will always return the last filename used for a load or save operation.

The third is similar but convert the file to *cvs* format while the fourth also specifies that the delimiter is a tab character.

2.5 Plotting Data

Data plotting and visualisation is handled by the **PlotFile** sub-class of **DataFile**. The purpose of the methods detailed here is to provide quick and convenient ways to plot data rather than providing publication ready figures.

```
>>> import Stoner.Plot as plot
>>> p=plot.PlotFile(d)
```

The first line imports the **Stoner.Plot** module. Strictly, this is unnecessary as the **Plot** module's namespace is imported when the **Stoner** package as a whole is imported. The second line creates an instance of the **PlotFile** class. **PlotFile** inherits the constructor method of **DataFile** and so all the variations detailed above work with **PlotFile**. In particular, the form shown in the second line is a easy way to convert a **DataFile** instance to a **PlotFile** instance for plotting.

2.5.1 Plotting 2D data

x-y plots are produced by the **PlotFile.plot_xy** method:

```
>>> p.plot_xy(column_x, column_y)
>>> p.plot_xy(column_x, [y1,y2])
>>> p.plot_xy(x,y,'ro')
>>> p.plot_xy(x,[y1,y2],['ro','b-'])
>>> p.plot_xy(x,y,title='My Plot')
>>> p.plot_xy(x,y,figure=2)
>>> p.plot_xy(x,y,plotter=pyplot.semilogy)
```

The examples above demonstrate several use cases of the **plot_xy** method. The first parameter is always the x column that contains the data, the second is the y-data either as a single column or list of columns. The third parameter is the style of the plot (lines, points, colours *etc*) and can either be a list if the y-column data is a list or a single string. Finally additional parameters can be given to specify a title and to control which figure is used for the plot. All matplotlib keyword parameters can be specified as additional keyword arguments and are passed through to the relevant plotting function. The final example illustrates a convenient way to produce log-linear and log-log plots. By default, **plotxy** uses the **pyplot.plot** function to produce linear scaler plots. There are a number of useful plotter functions that will work like this:

pyplot.semilogx, pyplot.semilogy These two plotting functions will produce log-linear plots, with **semilogx** making the x-axes the log one and **semilogy** the y-axis.

pyplot.loglog Like the semi-log plots, this will produce a log-log plot.

pyplot.errorbar this particularly useful plotting function will draw error bars. The values for the error bars are passed as keyword arguments, *xerr* or *yerr*. In standard matplotlib, these can be numpy arrays or constants. **PlotFile.plot_xy** extends this by intercepting these arguments and offering some short cuts:

```
>>> p.plot_xy(x,y,plotter=errorbar,yerr='dResistance',
              xerr=[5,'dTemp+'])
```

This is equivalent to doing something like:

```
>>> p.plot_xy(x,y,plotter=errorbar,
              yerr=p.column('dResistance'),
              xerr=[p.column(5),p.column('dTemp+')])
```

If you actually want to pass a constant to the *x/yerr* keywords you should use a float rather than an integer.

The X and Y axis label will be set from the column headers.

2.5.2 Plotting 3D Data

A number of the measurement rigs will produce data in the form of rows of *x, y, z* values. Often it is desirable to plot these on a surface plot or 3D plot. The **PlotFile.plot_xyz** method can be used for this.

```
>>> p.plot_xyz(col_x,col_y,col_z)
>>> p.plot_xyz(col_x,col_y,col_z,cmap=matplotlib.cm.jet)
>>> p.plot_xyz(col_x,col_y,col_z,plotter=pyplot.pcolor)
```

By default the **plot_xyz** will produce a 3D surface plot with the z-axis coded with a rainbow colourmap (specifically, the matplotlib provided *matplotlib.cm.jet* colourmap. This can be overridden with the *cmap* keyword parameter. If a simple 2D surface plot is required, then the *plotter* parameter should be set to a suitable function such as **pyplot.pcolor**.

Like **plot_xy**, a *figure* parameter can be used to control the figure being used and any additional keywords are passed through to the plotting function. The axes labels are set from the corresponding column labels.

Alternatively, if your data is already in the form of a matrix, you can use the **PlotFile.plot_matrix** method:

```
>>> p.plot_matrix()
>>> p.plot_matrix(xvals,yvals,rectang,title="Title",xlabel="X Axis",
                  ylabel="Y Axis",zlabel="Z Axis",cmap=matplotlib.cm.jet)
>>> p.plot_matrix(plotter=pyplot.pcolor,figure=False)
```

The first example just uses all the default values, in which case the matrix is assumed to run from the 2nd column in the file to the last and over all of the rows. The x values for each row are found from the contents of the first column, and the y values for each column are found from the column headers interpreted as a floating pint number. The colourmap defaults to the built in ‘jet’ theme. The x axis label is set to be the column header for the first column, the y axis label is set either from the meta data item “ylabel” or to “Y Data”. Likewise the z axis label is set from the corresponding metadata item or defaults to “Z Data”. In the second form these parameters are all set explicitly. The *xvals* parameter can be either a column index (integer or string) or a list, tuple or numpy array. The *yvals* parameter can be either an row number (integer) or list,tuple or numpy array. Other parameters (including *plotter*, *figure* etc) work as for the **PlotFile.plot_xyz** method. The *rectang* parameter is used to select only part of the data array to use as the matrix. It may be 2-tuple in which case it specifies just the origin as (row,column) or a 4-tuple in which case the third and forth elements are the number of rows and columns to include. If *xvals* or *yvals* specify particular column or rows then the origin of the matrix is moved to be one column further over and one row further down (*i.e.* the matrix is to the right and below the columns and rows used to generate the x and y data values). The final example illustrates how to generate a new 2D surface plot in a new window using default matrix setup.

2.5.3 Getting More Control on the Figure

It is useful to be able to get access to the matplotlib figure that is used for each **PlotFile** instance. The **PlotFile.fig** attribute can do this, thus allowing plots from multiple **PlotFile** instances to be combined in a single figure.

```
>>> p1.plot_xy(0,1,'r-')
>>> p2.plot_xy(0,1,'bo',figure=p1.fig)
```

Likewise the **PlotFile.axes** attribute returns the current axes object of the current figure in use by the **PlotFile** instance.

There’s a couple of extra methods that just pass through to the pyplot equivalents:

```
>>> p.draw()
>>> p.show()
```

2.6 Manipulating and Analysing Data

Curve fitting, data manipulation, and other analysis functions is handled by a sub-class of the DataFile object – AnalyseFile

```
>>> import Stoner.Analysis as Analysis
>>> a=Analysis.AnalyseFile('Data')
>>> a2=Analysis.AnalyseFile()
>>> a2=d
>>> a3=Analysis.AnalyseFile(d)
```

The first line imports the `AnaylseFile` class. Since the `AnalyseFile` is a child class of `DataFile`, everything you can do with a `DataFile` also works with an `AnalyseFile` object. The next two lines demonstrate creating a blank `AnalyseFile` and then copying all of the data, metadata and column headings from an existing `dataFile` object. The final variant shows how to cast one child-class of **DataFile** into another – in this case an **AnalyseFile**.

2.7 Manipulating Data

Several methods are provided to assist with common data fitting and preparation tasks, such as normalising columns, adding and subtracting columns.

```
>>> a.normalise('data','reference',header='Normalised Data',replace=True)
>>> a.normalise(0,1)
>>> a.normalise(0,3.141592654)
>>> a.normalise(0,a2.column(0))
```

The *normalise* method simply divides the data column by the reference column. By default the *normalise* method replaces the data column with the new (normalised) data and appends “(norm)” to the column header. The keyword arguments `header` and `replace` can override this behaviour. The third variant illustrates normalising to a constant (note, however, that if the second argument is an integer it is treated as a column index and not a constant). The final variant takes a 1D array with the same number of elements as rows and uses that to normalise to. A typical example might be to have some baseline scan that one is normalising to.

```
>>> a.subtract('A','B',header="A-B",replace=True)
>>> a.subtract(0,1)
>>> a.subtract(0,3.141592654)
>>> a.subtract(0,a2.column(0))
```

As one might expect from the name, the *subtract* method subtracts the second column from the first. Unlike *normalise* the first data column will not be replaced but a new column inserted and a new header (defaulting to column header 1 - column header 2) will be created. This can be overridden with the `header` and `replace` keyword arguments. The next two variants of the *subtract* method work in an analogous manner to the *normalise* methods. Finally the *add* method allows one to add two columns in a similar fashion:

```
>>> a.add('A','B',header='A plus B',replace=False)
>>> a.add(0,1)
>>> a.add(0,3.141592654)
>>> a.add(0,a2.column(0))
```

For completeness we also have:

```
>>> a.divide('A','B',header='A/B', replace=True)
>>> a.multiply('A','B',header='A*B', replace=True)
```

with variants that take either a 1D array of data or a constant instead of the B column index.

2.8 Curve Fitting

2.8.1 Simple polynomial Fits

Simple least squares fitting of polynomial functions is handled by the **AnalyseFile.polyfit** method:

```
>>> a.polyfit(column_x,column_y,polynomial_order, bounds=lambda x, y:True,
               result="New Column")
```

This is a simple pass through to the numpy routine of the same name. The *x* and *y* columns are specified in the first two arguments using the usual index rules for the Stoner package. The routine will fit multiple columns if *column_y* is a list or slice. The *polynomial_order* parameter should be a simple integer greater or equal to 1 to define the degree of polynomial to fit. The bounds function follows the same rules as the bounds function in **DataFile.search** to restrict the fitting to a limited range of rows. The method returns a list of co-efficients with the highest power first. If *column_y* was a list, then a 2D array of co-efficients is returned.

If *result* is specified then a new column with the header given by the result parameter will be created and the fitted polynomial evaluated at each point.

2.8.2 Simple function fitting

For more general curve fitting operations the **AnalyseFile.cruve_fit** method can be employed. Again, this is a pass through to the numpy routine of the same name.

```
>>> a.curve_fit(func, xcol, ycol, p0=None, sigma=None,
               bounds=lambda x, y: True, result="New column" )
```

The first parameter is the fitting function. This should have prototype *y=func(x,p[0],p[1],p[2]...)* where *p* is a list of fitting parameters. The *p0* parameter contains the initial guesses at the fitting parameters, the default value is 1. *xcol* and *ycol* are the *x* and *y* columns to fit. This method cannot handle multiple *y* columns. *sigma*, if present, provides the weightings for each datapoint and so should also be an array of the same length as the *x* and *y* data. Fianlly, the bounds function can be used to restrict the fitting to only a subset of the rows of data.

AnalyseFile.curve_fit returns a list of two arrays [*popt*,*pcov*] where *popt* is an array of the optimal fitting parameters and *pcov* is a 2D array of the co-variances between the parameters.

If *result* is specified then a new column with the header given by the result parameter will be created and the fitted function evaluated at each point.

2.8.3 Fitting with limits

For cases where one requires more flexibility in fitting data, in particular where the fitting parameters are constrained, the **AnalyzeFile.mpfir** method is provided. This is a pass through to the **mpfir** module.

```
>>> a.mpfir(func, xcol, ycol, p_info, func_args=dict(), sigma=None,
            bounds=lambda x, y: True, **mpfir_kargs )
```

In this case, the *func* argument takes a slightly different prototype: **def func(x,parameters, **func_args)** where *parameters* is a list of the fitting parameters and *func_args* provides a dictionary of fixed *i.e.* non-fitting parameters. *xcol*

and *ycol* are the column indices for the x and y data, *bounds* is a bounding function to select only those rows to use for fitting the function, and *sigma* are the weightings for each datapoint. The remaining arguments are a dictionary of keywords to pass through to the **mpfit** routine and *p_info* which is a list of dictionaries which is used to control the parameters in the fit. This described below.

p_info contains one element for each parameter used to fit the data. Each element is a dictionary with the following keys:

value the starting parameter value (but see the START_PARAMS parameter for more information).

fixed a boolean value, whether the parameter is to be held fixed or not. Fixed parameters are not varied by MPFIT, but are passed on to MYFUNCT for evaluation.

limited a two-element boolean array. If the first/second element is set, then the parameter is bounded on the lower/upper side. A parameter can be bounded on both sides. Both LIMITED and LIMITS must be given together.

limits a two-element float array. Gives the parameter limits on the lower and upper sides, respectively. Zero, one or two of these values can be set, depending on the values of LIMITED. Both LIMITED and LIMITS must be given together.

parname a string, giving the name of the parameter. The fitting code of MPFIT does not use this tag in any way. However, the default iterfunct will print the parameter name if available.

step the step size to be used in calculating the numerical derivatives. If set to zero, then the step size is computed automatically. Ignored when AUTODERIVATIVE=0.

mpside the sidedness of the finite difference when computing numerical derivatives. This field can take four values:

0 one-sided derivative computed automatically

1 one-sided derivative $(f(x+h) - f(x))/h$

-1 one-sided derivative $(f(x) - f(x-h))/h$

2 two-sided derivative $(f(x+h) - f(x-h))/(2 * h)$

Where H is the STEP parameter described above. The "automatic" one-sided derivative method will chose a direction for the finite difference which does not violate any constraints. The other methods do not perform this check. The two-sided method is in principle more precise, but requires twice as many function evaluations. Default: 0.

mpmaxstep the maximum change to be made in the parameter value. During the fitting process, the parameter will never be changed by more than this value in one iteration.

A value of 0 indicates no maximum. Default: 0.

tied a string expression which "ties" the parameter to other free or fixed parameters. Any expression involving constants and the parameter array P are permitted. Example: if parameter 2 is always to be twice parameter 1 then use the following: `parinfo(2).tied = '2 * p(1)'`. Since they are totally constrained, tied parameters are considered to be fixed; no errors are computed for them. [NOTE: the PARNAME can't be used in expressions.]

mpprint if set to 1, then the default iterfunc will print the parameter value. If set to 0, the parameter value will not be printed. This tag can be used to selectively print only a few parameter values out of many.
 Default: 1 (all parameters printed)

2.9 More AnalyseFile Functions

2.9.1 Basic Data Inspection

```
>>> a.max(column)
>>> a.min(column)
>>> a.max(column,bounds=lambda x,y:y[2]>1 and y[2]<10)
>>> a.min(column,bounds=lambda x,y:y[2]>1 and y[2]<10)
```

Hopefully all of the above are fairly obvious ! In the last two cases, one can use a function to limit the search to particular rows (*e.g.* to search for the maximum y value subject to some constraint in x). One important point to note is that the routines return a tuple of two numbers, the maximum (or minimum) and the row number where the maximum or minimum was found.

There are a couple of related functions to help here:

```
>>> a.span(column)
>>> a.span(column, bounds=lambda x,y:y[2]>100)
>>> a.clip(column,(max_v,min_v)
>>> a.clip(column,b.span(column))
```

The **span** method simply returns a tuple of minimum and maximum values within either the whole column or bounded data. Internally this is just calling the **max** and **min** methods. The **clip** method deletes rows for which the specified column as a value that is either larger or smaller than the maximum or minimum value within the second argument. This allows one to specify either a tuple – *e.g.* the result of the **span** method, or a complete list as in the last example above. Specifying a single float would have the effect of removing all rows where the column didn't equal the float value. This is probably not a good idea...

It is worth pointing out that these functions will respect the existing mask on the data unless the bounds parameter is set, in which case the mask is temporarily discarded in favour of one generated from the bounds expression. This can be worked around, however, as the parameter passed to the bounds function is itself a masked array and thus one can include a test of the mask in the bounds function:

```
>>> a.span(column,bounds=lambda x,y:y[2]>10 or not numpy.any(y.mask))
```

2.9.2 Thresholding and Interpolating Data

```
>>> a.threshold(col, threshold, rising=True, falling=False)

>>> a.interpolate(newX,kind='linear' )
```

2.9.3 Smoothing and Differentiating Data

2.9.4 Peak Finding

2.10 Working with Lots of Files

A common case is that you have measured lots of data curves and now have a large stack of data files sitting in a tree of folders on disc and now need to process all of them with some code. The **DataFolder** class is designed to make it easier to process lots of files.

2.10.1 Getting a List of Files

The first thing you probably want to do is to get a list of data files in a directory (possibly including its subdirectories) and probably matching some sort of filename pattern.

```
>>> from Stoner.Folders import DataFolder
>>> f=DataFolder(pattern='*.dat')
```

In this very simple example, the **DataFolder** class is imported in the first line and then a new instance *f* is created. The optional *pattern* keyword is used to only collect the files with a .dat extension. In this example, it is assumed that the files are readable by **DataFile**, if they are in some other format then the *type* keyword can be used:

```
>>> from Stoner.FileFormats import XRDFFile
>>> f=DataFolder(type=XRDFFile,pattern='*.dql')
```

To specify a particular directory to look in, simply give the directory as the first argument - otherwise the current directory will be used.

```
>>> f=DataFolder('/home/phygbu/Data',pattern='*.tdi')
```

By default the **DataFolder** constructor will perform a recursive directory listing of the working folder. The resulting list of files can be accessed via the *files* attribute:

```
>>> f.files
```

If you don't want the file listing to be recursive, this can be suppressed by using the *recursive* keyword argument and the file listing can be suppressed altogether with the *nolist* keyword:

```
>>> f=DataFolder(pattern='*.dat',recursive=False)
>>> f2=DataFolder(nolist=True)
```

The current root directory and pattern are stored in the *directory* and *pattern* keywords and the **getlist** method can be used to force a new listing of files.

```
>>> f.directory='/home/phygbu/Data'
>>> f.pattern='*.txt'
>>> f.getlist()
```

Sometimes a more complex filename matching mechanism than simple “globbing” is useful. The *pattern* keyword can also be a compiled regular expression:

```
>>> import re
>>> p=re.compile('i10-\d*.dat')
>>> f=DataFolder(pattern=p)
>>> p2=re.compile('i10-(?P<run>\d*)')
>>> f=DataFolder(pattern=p)
>>> f[0]['run']
```

The second case illustrates a useful feature of regular expressions - they can be used to capture parts of the matched pattern – and in the python version, one can name the capturing groups. In both cases above the **DataFolder** has the same file members (basically these would be runs produced by the i10 beamline at Diamond), but in the second case the run number (which comes after “i10-” would be captured and presented as the “run” parameter in the metadata when the file was read. **Note that the files**



are not modified - the extra metadata is only added as the file is read by the DataFolder . The loading process will also add the metadata key “Loaded From” to

the file which will give you a note of the filename used to read the data.

Finally, akin to **DataFile** you can force a dialog box to select a directory by passing *False* into the constructor or getlist methods in place of a directory name.

2.10.2 Doing Something With Each File

A **DataFolder** is an object that you can iterate over, loading the **DataFile** type object for each of the files in turn. This provides an easy way to run through a set of files, performing the same operation on each:

```
>>> folder=DataFolder(pattern='*.tdi')
>>> for f in folder:
....     f=AnalyseFile(f)
....     f.normalise('mac116','mac119')
....     f.save()
```

or even more compacts:

```
>>>[f.normalise('mac116','mac119').save() for f in
      DataFolder(pattern='*.tdi',type=AnalyseFile)]
```

DataFolder is also indexable and has a length:

```
>>> f=DataFolder()
>>> len(f)
>>> f[0]
```

2.10.3 Sorting, Filtering and Grouping Data Files

The order of the files in a **DataFolder** is arbitrary. If it is important to process them in a given order then the *sort* method can be used:

```
>>> f.sort()
>>> f.sort('temperature')
>>> f.sort('Temperature',reverse=True)
>>> f.sort(lambda x:len(x))
```

The first variant simply sorts the files by filename. The second and third variants both look at the “temperature” metadata in each file and use that as the sort key. In the third variant, the *revers* keyword is used to reverse the order of the sort. In the final variant, each file is loaded in turn and the supplied function is called and evaluated to find a sort key.

The **filter** method can be used to prune the list of files to be used by the **DataFoler**:

```
>>> f.filter('[ab]*.dat')
>>> import re
>>> f.filter(re.compile('i10-\\d*\\.dat'))
>>> f.filter(lambda x: x['Temperature']>150)
>>> f.filter(lambda x: x['Temperature']>150, invert=True)
```

The first form performs the filter on the filenames (using the standard python fnmatch module). One can also use a regular expression as illustrated in the second example – although unlike using the *pattern* keyword in **getlist**, there is no option to capture metadata (although one could then subsequently set the pattern to achieve this). The third variant calls the supplied function, passing the current file as a **DataFile** object in each time. If the function evaluates to be True then the file is kept. The *invert* keyword is used to invert the sense of the filter (a particularly silly example here, since the greater than sign could simply be replaced with a less than or equals sign !).

One of the more common tasks is to group a long list of data files into separate groups according to some logical test – for example gathering files with magnetic field sweeps in a positive direction together and those with magnetic field in a negative direction together. The **group** method provides a powerful way to do this. Suppose we have a series of data curves taken at a variety of temperatures and with three different magnetic fields:

```
>>> f.group{'temperature'}
>>> f.group(lambda x:"positive" if x['B-Field']>0 else "negative")
>>> f.group(['temperature', lambda x:"positive"
            if x['B-Field']>0 else "negative"])
>>> f.groups
```

The **group** method splits the files in the **DataFolder** into several groups each of which share a common value of the argument supplied to the **group** method. A group is itself another instance of the **DataFolder** class. Each **DataFolder** object maintains a dictionary called *groups* whose keys are the distinct values of the argument of the **group** methods and whose values are **DataFolder** objects. So, if our **DataFolder** *f* contained files measured at 4.2, 77 and 300K and at fields of 1T and -1T then the first variant would create 3 groups: 4.2, 77 and 300 each one of which would be a **DataFolder** object containing the files measured at those temperatures. The second variant would produce 2 groups – “positive” containing the files measured with magnetic field of 1T and “negative” containing the files measured at -1T. The third variant then goes one stage further and would produce 3 groups, each of which in turn had 2 groups. The groups are accessed via the *group* attribute:

```
>>> f.groups[4.2].groups["positive"].files
```

would return a list of the files measured at 4.2K and 1T.

In some cases you will want to work with sets of files coming from different groups in order. For example, if above we had a sequence of 10 data files for each field and

temperature and we wanted to process the positive and negative field curves together for a given temperature in turn. In this case the `zip_groups` method can be useful.

```
>>> f.groups[4.2].zip_groups(['positive','negative'])
```

This would return a list of tuples of `DataFile` objects where the tuples would be the first positive and first negative field files, then the second of each, then third of each and so. This presupposes that the files started of sorted by some suitable parameter (*e.g.* a gate voltage).

3 Cookbook

This section gives some short examples to give an idea of things that can be done with the Stoner python module in just a few lines.

3.1 Extract X-Y(Z) from X-Y-Z data

In a number of measurement systems the data is returned as 3 parameters X, Y and Z and one wishes to extract X-Y as a function of constant Z. For example, $I - V$ sweeps as a function of gate voltage V_G . Assuming we have a data file with columns *Current*, *Voltage*, *Gate*:

```
>>> d=DataFile('data.txt')
>>> t=d
>>> for gate in d.unique('Gate'):
>>>     t.data=d.search('Gate',gate)
>>>     t.save('Data Gate='+str(gate)+'.txt')
```

The first line opens the data file containing the $I - V(V_G)$ data. The second creates a temporary copy of the `DataFile` object - ensuring that we get a copy of all metadata and column headers. The **for** loop iterates over all unique values of the data in the gate column and then inside the for loop, searches for the corresponding $I - V$ data, sets it as the data of the temporary `DataFile` and then saves it.

3.2 Mapping X-Y-Z data to Z(X,Y) data

In a similar fashion to the previous section, where data has been recorded with fixed values of X and Y *e.g.* I measured for fixed V and V_G , it can be useful to map the data to a matrix.

```
>>> d=DataFile('Data,.txt')
>>> t=d
>>> for gate in d.unique('Gate'):
>>>     t=t+d.search('Gate',gate)[: ,d.find_col('Current')]
>>> t.column_headers=['Bias='+str(x) for x in d.unique('Voltage')]
>>> t.add_column(d.unique('Gate'),'Gate Voltage',0)
```

The start of the script follows the previous section, however this time in the for loop the addition operator is used to add a single row to the temporary `DataFile` *t*. In this case we are using the utility method **DataFile.find_col** to find the index of the column with the current data. After the for loop we set the column headers in *t* and then insert an additional column at the start with the gate voltage values.

The matrix generated by this code is suitable for feeding directly into **PlotFile.plot_matrix()**, however, the same plot could be generated directly from the **PlotFile.plot_xyz()** method too.

4 Developer's Guide

This section provides some notes and guidance on extending the Stoner Package.

4.1 Adding New Data File Types

The first question to ask is whether the data file format that you are working with is one that others in the group will be interested in using. If so, then the best thing would be to include it in the **fileFormats** module in the package, otherwise you should just write the class in your own script files. In either case, develop the class in your own script files first.

The best way to implement handling a new data format is to write a new subclass of **DataFile**:

```
class NewInstrumentFile(DataFile):
    """Extends DataFile to load files from somewhere else

    Written by Gavin Burnell 11/3/2012"""
```

A document string should be provided that will help the user identify the function of the new class (and avoid using names that might be commonly replicated !). Only one method needs to be implemented: a new *load* method. The *load* method should have the following structure:

```
def load(self,filename=None,*args):
    """Just call the parent class but with the right parameters set"""
    if filename is None or not filename:
        self.get_filename('r')
    else:
        self.filename = filename
```

then follows the code to actually read the file. It must at the very least provide a column header for every column of data and read in as much numeric data as possible and it should read as much of the meta data as possible. The function terminates by returning a copy of the current object:

```
    return self
```

One useful function for reading metadata from files is *self.metadata.string_to_type()* which will try to convert a string representation of data into a sensible Python type.

There is one global attribute that can be used to tweak the automatic file importing code.

```
f.priority=32
```

When the subclasses are tried to see if they can load an undetermined file, they are tried in order of priority. If your load code can make a positive determination that it has the correct file (*e.g.* by looking for some magic combination of characters at the start of the file) and can throw an exception if it tries loading an incorrect file, then you can

give it a lower priority number to force it to run earlier. Conversely if your only way of identifying your own files is seeing they make sense when you try to load them and that you might partially succeed with a file from another system (as can happen if you have a tab or comma separated text file), then you should raise the priority number. Currently **DataFolder** defaults to 32, **CSVFile** and **BigBlueFile** have values of 128 and 64 respectively.

If you need to write any additional methods please make sure that they have DoxyGen document strings so that the API documentation is picked up correctly.